

**URGENT FIELD SAFETY NOTICE**

No. 1226348-9/10/14-001R

September 12, 2014

**Voluntary Recall Notification:  
EDS 3™ CSF External Drainage System****PLEASE DISTRIBUTE THIS INFORMATION TO  
CLINICIANS WHO USE THIS DEVICE**

Dear OR Manager, Materials Manager, ICU/Neuro ICU Manager:

Please be aware that Codman Neuro is initiating a voluntary recall of all lots of EDS 3 CSF External Drainage Systems, with expiration on or before August 2017 (2017-08, please see attachment). The system is indicated for draining cerebrospinal fluid (CSF) and other fluids of similar characteristics as a means of reducing intracranial pressure and CSF volume when the insertion of a permanent, internal shunt is not indicated. Please be aware that products in your possession may be affected.

**Affected Product (all lots)**

<b>Code</b>	<b>Description</b>
82-1730	EDS 3CSF External Drainage System with Ventricular Catheter
82-1731	EDS 3CSF External Drainage System (no Ventricular Catheter)
82-1738	Lumbar Drainage Catheter Kit II with EDS 3System

This recall is being initiated because the tubing within the system that drains CSF may leak or disconnect from the joints. Leakage and tube separations may result in over- or under-drainage of CSF from the ventricular system or introduction of air into the ventricular system (pneumocephalus). This may result in collapsed ventricles, subdural bleeding, or an inability to properly control elevated intracranial pressure. The tubing disconnection or leakage may also increase the risk of ventriculitis. If undetected or untreated each of these events may cause severe brain injury, which may lead to coma, stroke or death.

**Actions for Clinicians Treating Patients with the EDS 3 System**

These systems are most often used on neurocritical care floors and these issues are likely to be detected immediately. For patients currently being managed with the EDS 3 System, the system

should be replaced immediately. In the case where no substitute drainage system is immediately available, the EDS 3 System may continue to be used until an alternative product can be obtained. Manipulation of tubing should be minimized and extra vigilance is required for early detection of leakage and/or disconnection.

As of September 9, 2014, approximately 0.051 percent of complaints reported a failed connection from the joint and 0.030 percent reported leakage according to the most recent complaint data, but some customers' complaint rates have been as high as 2.6 percent.

**The relevant National Competent Authorities are informed of this action. Please provide this notice to any neurosurgeons or other clinicians at your facility who are managing patients with the EDS 3 System, and to other organizations to which products may have been transferred.**

We request that you immediately check all inventory to determine if you have affected product. Please use the included instructions to report your inventory status and return affected product.

We apologize for the inconvenience this recall may cause. At Codman Neuro, we are dedicated to delivering products that meet the highest quality standards. We regret the need to undertake this voluntary recall and we are committed to resolving this issue.

If you have any questions or concerns regarding this notification, please contact your local Codman Neuro Representative or Scientific and Medical Affairs at [SciMedAffairs@its.jnj.com](mailto:SciMedAffairs@its.jnj.com) or +1 (866) 685-7325. Thank you for your cooperation.

Sincerely,



J. Thomas Megerian MD, PhD  
Vice President - Strategic Medical Affairs and Medical Sciences

**CODMAN NEURO**



**FIELD SAFETY NOTICE No. 1226348-9/10/14-001R  
INSTRUCTIONS FOR REPORTING INVENTORY STATUS, RETURNING PRODUCT,  
AND OBTAINING CREDIT**

You are receiving this voluntary recall notification because our records indicate that you are the recipient of one or more EDS 3 systems affected by this recall. Please follow the instructions in this field safety notice for the products identified in this recall. Please complete the enclosed Acknowledgement Form and fax the completed form to: **001-508-977-6665**. **It is important that we receive this acknowledgement form, even if you have no affected product in your inventory.**

You may also e-mail the form to us at [Codmanproductcomplaints@dpyus.jnj.com](mailto:Codmanproductcomplaints@dpyus.jnj.com)

**PRODUCT RETURNS:**

All affected products are to be returned to the address below for credit.

*GMED Healthcare  
EDC Quality Dept  
Rue de Luxembourg 5  
ZI Trazegnies  
BE - 6180 Courcelles  
Belgium*

**TEL: 32-7-146-9404**

Once again, it is very important that you complete the acknowledgement form, as soon as possible, even if you do not have any affected product inventory in your possession. You may fax the completed acknowledgement form to the U.S. at **001-508-977-6665**. You may also e-mail the form to us at [Codmanproductcomplaints@dpyus.jnj.com](mailto:Codmanproductcomplaints@dpyus.jnj.com). If you have any questions or concerns in regards to this recall, please contact your Codman Neuro consultant.

# CODMAN NEURO



## Field Safety Notice Acknowledgement Form FSN: 1226348-9/10/14-001R

This signed form acknowledges receipt of the recall letter issued by Codman Neuro regarding the EDS 3™ CSF External Drainage System. Not completing this form in its entirety can cause a delay in processing your credit. Please return this completed form to your local Codman Neuro Consultant or fax to: < NOTE: Add Local Fax Number Here >

We have checked our current inventory for affected products listed in the recall announcement and indicated our inventory levels below:

- We acknowledge the receipt of this information but **do not have** the affected EDS 3 CSF External Drainage System products in stock.
- We acknowledge the receipt of this information and the products listed below **are in stock** and will be returned to Codman Neuro for dollar credit.
- We have affected product, but it is currently in use with no available alternative, so we are retaining product to address emergent patient need. We will return affected product once an alternative product is available.

<i>Product Code</i>	<i>Lot Number</i>	<i>Quantity</i>

Hospital Name/Location: \_\_\_\_\_

Print Name: \_\_\_\_\_

Phone Number: \_\_\_\_\_

---

**Authorized Signature/Date**

# CODMAN NEURO



## ATTACHMENT 1

### EDS 3™ CSF External Drainage System (Sample) Product Label

Codman & Shurtleff, Inc  
325 Paramount Drive  
Raynham, MA 02767-0350, USA

EC REP Johnson & Johnson Medical, Ltd  
Pinewood Campus, Nine Mile Ride  
Wokingham, RG40 3EW, UK

**MADE IN** Switzerland of US and foreign components

For patent information about this product, go to  
[www.depu.com/patentmarking](http://www.depu.com/patentmarking)



(01)10886704040903(20)99



(17)991231(10)SAMPLE




LCN2 REV. H

**Expiration date information located here. All product affected expires on or before 2017-08.**

**REF** 82-1730  YYYY-MM  YYYY-MM

**LOT** SAMPLE **QTY** 1

**CODMAN® EDS 3™**  
CSF External Drainage System with Ventricular Catheter  
(REF 82-1735)



Système de drainage externe de LCR avec cathéter ventriculaire  
(REF 82-1735)

Externes Liquordrainagesystem mit Ventrikelkatheter (REF 82-1735)

Extern liquor-drainagesysteem met ventrikelkatheter (REF 82-1735)

Sistema di drenaggio esterno per liquido cerebrospinale con  
catetere ventricolare (REF 82-1735)

Sistema de drenaje externo de LCR con catéter ventricular (REF  
82-1735)

Sistema de drenagem externa de LCR com cateter ventricular (REF.  
82-1735)

Eksternt cerebrospinalvæske-drænagesystem med ventrikulært  
kateter (REF 82-1735)

Externt likvordrænagesystem med ventrikelkateter (REF 82-1735)

Ulkoinen CSF-aivo-selkäydinnesteen dreneerausjärjestelmä ja  
aivokammiokatetri (REF 82-1735)

Σύστημα εξωτερικής παροχέτευσης ΕΝΥ με κοιλιακό καθετήρα (REF  
82-1735)

A cerebrospinális folyadékot elvezető, ventricularis katéterrel  
felszerelt külső drenázsrendszer (REF 82-1735)