

Report Form Manufacturer's Field Safety Corrective Action Report

Medical Devices Vigilance System (MEDDEV 2.12/1 rev 8)

v.01.13

1. Administrative information	
To which NCA(s) is this report being sent? Agency for Medicines and Medical Devices of Bosnia and Herzegovina	
Type of report <input checked="" type="checkbox"/> Initial report <input type="checkbox"/> Follow up report <input type="checkbox"/> Final report	
Date of this report 8 January 2014	
Reference number assigned by the manufacturer FA599	
FSCA reference number assigned by NCA	
Incidence reference number assigned by NCA	
Name of the co-ordinating national competent authority (if applicable)	
2. Information on submitter of the report	
Status of submitter <input type="checkbox"/> Manufacturer <input checked="" type="checkbox"/> Authorised representative within EEA, Switzerland and Turkey <input type="checkbox"/> Others (identify the role): Manufacturer legal entity in Switzerland	
3 Manufacturer information	
Name Medtronic Inc.	
Contact name	
Address 710 Medtronic Parkway	
Postcode MN 55432-5604	City Minneapolis
Phone (+1-763) 514-4000	Fax (+1-763) 514-4879
E-mail	Country USA
4 Authorised representative information	
Name Medtronic BV	
Contact name Jean-Charles Moreau	
Address Earl Bakkenstraat 10	

Postcode 6422PJ	City Heerlen
Phone +33 1 55 38 60 35	Fax +33 1 55 38 18 00
E-mail vigilance.eu@medtronic.com	Country NL
5 National contact point information	
National contact point name	
Name of the contact person	
Address	
Postal code	City
Phone	Fax
E-mail	Country
6 Medical device information	
Class	
<input type="checkbox"/> AIMD Active implants	<input type="checkbox"/> IVD Annex II List A
<input checked="" type="checkbox"/> MDD Class III	<input type="checkbox"/> IVD Annex II List B
<input type="checkbox"/> MDD Class IIb	<input type="checkbox"/> IVD Devices for self-testing
<input type="checkbox"/> MDD Class IIa	<input type="checkbox"/> IVD General
<input type="checkbox"/> MDD Class I	
Nomenclature system (preferable GMDN) GMDN	Nomenclature code 15870
Nomenclature text Prosthesis, cardiac valve, biological	
Commercial name/brand name/make Medtronic Mosaic™ Aortic Bioprosthesis	
Model number 305	Catalogue number
Serial number(s) All	lot/batch number(s)
Device Manufacturing date NA	Expiry date NA
Software version number (if applicable)	
Accessories/associated device (if applicable) n/a	
Notified body (NB) ID- number 0123	
7 Description of FSCA	
<p>Device description:</p> <p>The Mosaic Aortic Bioprosthesis is indicated for the replacement of malfunctioning native or prosthetic aortic valves and is available in sizes 19mm through 29mm.</p> <p>Description of Issues:</p> <p>There have been a total of 679 reported events for high gradient (inclusive of reports of high gradient and stenosis) from January 2000 through October 2013. The complaint trend has been stable for the last five years, with approximately 80 reported events each year, with the exception from May 2010 to April 2011, which had 116 reported events. High gradients and stenosis are considered chronic events, and as such are not expected to occur</p>	

within short implant durations without assignable physiologic and/or anatomical cause.

There have been a total of 399 Mosaic valves explanted for high gradient and stenosis from May 2000 to October 2013 (349 Aortic/ 50 Mitral). Of the 349 Aortic valves explanted, 284 (Aortic) were returned to Medtronic for analysis. Based on total sales for Mosaic Aortic (224,427) this represents an explant rate of 0.16%. For the same period of time, for Mosaic Mitral™ (Model 310) the explant rate is 0.06% (50/83.105), which is in-line with the predicted rate.

Current root cause concludes that:

- The flexible stent of the Mosaic aortic valve may be susceptible to annular deformation if oversized with respect to the native aortic annulus (stent deformation).
- Deformation of the stent annulus can result in a high measured gradient. Echo measurements tend to overstate gradients compared with pressure transducer or catheter measurements. Overstated gradients tend to be more pronounced in valves exhibiting deformation. This phenomenon was clearly demonstrated on the bench using echo measurements (annular deformation can cause a higher gradient).
- In the absence of stent deformation, oversizing (upsizing) the valve can obstruct the flow area of the valve and/or cause altered leaflet kinematics resulting in an elevated gradient (bioprosthetic inflow obstruction).

The higher-than-expected transvalvular gradients (> 25 mmHg) of the Mosaic Aortic Bioprosthesis have occurred at a rate of 0.33% (3.3 reports per 1000 aortic implants). A subset of these valves was explanted within five years of implantation at an overall occurrence rate of 0.1% (1 explant per 1000 aortic implants).

The complaint analysis demonstrates that the current rates, based on cumulative data, and explants returned for analysis, are within the predicted risk, with the exception of Structural Valve Dysfunction resulting in Stenosis and Tissue Overgrowth resulting in Stenosis. The predicted risk for each of the failure modes (hazards) that can lead to the patient effect (or harm) of high gradient and/or stenosis leading to explant are "Broadly Acceptable" or "As Low as Reasonable Practical" (ALARP). The clinical performance of the valve based on the most recent Clinical Evidence Report, which includes a review of current literature from the last five years (Jan. 2007 to Mar. 2012) does not indicate or suggest that the overall safety and performance of the valve has changed relative to the predicted risk.

Description and justification of the action (corrective/preventive)

1. As a safety precaution measure, Medtronic is informing the customers (Physicians, Hospital Administrators, OR Managers and Risk Managers) of these issues.
2. Medtronic has modified the iEOA Mosaic Aortic Bioprosthesis sizing chart.
3. Medtronic has modified the current Mosaic Obturators/Sizers to align with the updated sizing chart. The IFU for the modified Mosaic Obturators/Sizers has been updated to include the new sizing chart (iEOA chart).
4. Medtronic will discontinue distribution of the previous Mosaic Obturators/Sizers and iEOA sizing chart.

Advice on actions to be taken by the distributor and the user:

Customers are recommended to dispose of existing iEOA sizing charts and begin to use the new iEOA sizing chart and the Mosaic Aortic Obturator/Sizer set.

Progress of FSCA, together with reconciliation data (Mandatory for a Final FSCA)

Attached please find

- Field Safety Notice (FSN) in English
 FSN in national language
 Others (please specify): Customer List in your country

FSN Status

- Draft
 Final

<p>Time schedule for the implementation of the different actions</p> <p>This FSCA is planned to be completed by 9 April 2014</p>
<p>These countries within the EEA and Switzerland and Turkey are affected by this FSCA</p> <p>- within the EEA, Switzerland and Turkey:</p> <p> <input checked="" type="checkbox"/> AT <input checked="" type="checkbox"/> BE <input checked="" type="checkbox"/> BU <input checked="" type="checkbox"/> CH <input checked="" type="checkbox"/> CY <input checked="" type="checkbox"/> CZ <input checked="" type="checkbox"/> DE <input checked="" type="checkbox"/> DK <input type="checkbox"/> EE <input checked="" type="checkbox"/> ES <input checked="" type="checkbox"/> FI <input checked="" type="checkbox"/> FR <input checked="" type="checkbox"/> GB <input checked="" type="checkbox"/> GR <input checked="" type="checkbox"/> HU <input checked="" type="checkbox"/> IE <input type="checkbox"/> IS <input checked="" type="checkbox"/> IT <input type="checkbox"/> LI <input checked="" type="checkbox"/> LT <input type="checkbox"/> LU <input checked="" type="checkbox"/> LV <input type="checkbox"/> MT <input checked="" type="checkbox"/> NL <input checked="" type="checkbox"/> NO <input checked="" type="checkbox"/> PL <input checked="" type="checkbox"/> PT <input type="checkbox"/> RO <input checked="" type="checkbox"/> SE <input checked="" type="checkbox"/> SI <input checked="" type="checkbox"/> SK <input checked="" type="checkbox"/> TR </p> <p>- Candidate Countries:</p> <p><input checked="" type="checkbox"/> HR</p> <p><input type="checkbox"/> All EEA, Candidate Countries, Switzerland and Turkey</p> <p>- Others:</p>
<p>8 Comments</p> <p>In Bosnia and Herzegovina customers are affected by this FSCA.</p>

I affirm that the information given above is correct to the best of my knowledge.

.....
Signature

Jean-Charles Moreau
Name

Heerlen
City

8 January 2014
Date

Submission of this report does not, in itself, represent a conclusion by the manufacturer and/or authorized representative or the national competent authority that the content of this report is complete or accurate, that the medical device(s) listed failed in any manner and/or that the medical device(s) caused or contributed to the alleged death or deterioration in the state of the health of any person.